Instructions for Use

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Bruksanvisning

INTENDED USE

The Altius Classic Central Venous Catheter Sets are sterile, single use devices containing a unique Altius central venous catheter and the required accessories for performing central venous catheterisation. Multiple variations of the catheter and the accessories are available to accommodate different patient anatomical requirements.

Central venous catheters are indicated for use in attaining access to the central venous system, and are used for the purpose of:

- administration of vasoactive/inotropic drugs,
- administration of incompatible medications,
- administration of hypertonic solutions including total parental nutrition,
- frequent blood sampling,
- blood transfusion,
- measurement of central venous pressure (CVP).

**WARNING**

- This product is sterile and for single use only. Do not re-use, reprocess or re-sterilise. Do not use catheter or accessories if any sign of damage is visible.
- Reprocessing or re-sterilisation of the set may damage the catheter and affect its integrity which may, when re-used, lead to severe deterioration in health and safety of patients.
- The catheter does not have any metallic components and can be exposed to various environmental conditions including thermal ignition source (during MRI) as long as no metal component is attached to it.
- Do not use catheter for contrast media injection otherwise the catheter may rupture and leak.
- The catheters should not be cut to alter their length.
- The catheters should be in place no longer than 30 days from the date of insertion.
- The insertion technique has a significant influence on the complications and outcome of the patient. Insertion must be performed by a competent and experienced catheter insertion team. Inexperienced personnel should not be permitted to perform the insertion except under the direct supervision of an experienced physician or surgeon.
- Be sure that you are familiar with the possible complications and emergency measures are known and available if occur.
- Patient requiring ventilation support is at greater risk of pneumothorax during sub-clavian vein cannulation.
- Do not advance the guidewire or catheter if unusual elastic resistance is encountered. Do not insert or withdraw the guidewire forcibly from any component. The wire could break or unravel, in which case both the catheter and guidewire must be removed simultaneously.
- In the rare event that a hub or connector separates from any component during insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism and remove the catheter immediately.

DEVICE DESCRIPTION

Altius central venous catheters are made from radiopaque polyurethane. The catheter tube can be single or multi-lumen. The catheter tip is tapered with softer material.

Kit components

Each Altius Classic Central Venous Catheter Set consists of a unique Altius central venous catheter and the following accessories:

- 1 × Introducer needle
- 1 × 5 ml guiding syringe or 1 × 3 ml standard syringe
- 1 × Nitinol guidewire in a guidewire advancer with cap
- 1 × Hydrophilic coated dilator
- 1 × Scalpel
- 1 × Fixation wing

Depending on the size, length and type of the catheter, different variations of accessories are provided in the sets. Dimensions of the components are indicated on the product label.

How supplied

All contents are supplied sterile and non-pyrogenic in unopened and undamaged packaging. The device is sterilised by ethylene oxide. Do not use catheter if sterile package has been damaged or has been opened. Do not re-use or re-sterilise the device.

Storage

Protect from direct sunlight, any source of heat, and moisture. Store and expose only within the specified temperature limit. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to expiration date on package label.

RELATIVE CONTRAINDICATIONS

The inserter of the catheter shall apply its clinical judgement prior to the insertion of the catheter. Instances where catheter might not be suitable, and not exhausted list, are:

- In patients with disorders
- When the presence of another device related infection, bacteraemia or septicaemia is known or suspected
- If severe chronic obstructive lung disease exists
- If previous episodes of venous thrombosis or vascular surgical procedure at the prospective insertion site have occurred
- Local tissue factors which may prevent proper devices stabilisation and/or access, like allergic reaction, or any dermatological disease

WHICH VEIN TO CANNULATE?

Possible insertion sites for central venous catheters are the right internal and external jugular veins, the left internal and external jugular veins, the subclavian veins and the femoral veins.

Table 1 – Patient evaluation prior to access placement

Consideration	Relevance
History of previous CVC	Previous placement of a CVC is associated with central venous stenosis.
History of pacemaker use	There is a correlation between pacemaker use and central venous stenosis.
History of diabetes mellitus	Diabetes mellitus is associated with damage to vasculature necessary for internal accesses.
History of anticoagulant therapy or any coagulation disorder	Abnormal coagulation may cause clotting or problems with haemostasis of access site.
History of vascular access	Previously failed vascular accesses will limit available sites for access; the cause of a previous failure may influence planned access if the cause is still present.

Consideration	Relevance
History of previous arm, neck or chest surgery/trauma	Vascular damage associated with previous surgery or trauma may limit viable access sites.

METHOD OF INSERTION

General preparation to obtain central access

The basic preparation and equipment which are required for venous cannulation are the same regardless of the route or technique chosen. Clinicians who insert central catheters should be taught the technique by an experienced colleague.

PRECAUTION

- Ultrasound should be used in the placement of catheters.
- The position of the tip of any central venous catheter should be verified by radiological means (e.g. X-ray) and monitored routinely as per institution policy.
- Do not use absolute alcohol or acetone based product on the catheter. 2% chlorhexidine or iodine based solution is recommended as antiseptic solution.
- It is not recommended to use antimicrobial ointments or solutions on catheters as it may cause degradation of catheter materials.
- Do not use sharp instruments near the extension line or tubing. Do not use scissors to remove dressing, as this could possibly cut or damage the catheter. Do not suture through any part of the catheter. Catheter tubing can tear when subjected to excessive force or rough edges.
- It is recommended that only luer lock (threaded) connections be used with the catheter (including syringes, bloodlines, IV tubing and injection caps).
- Repeated over-tightening of bloodlines, syringes and caps will reduce connector life and could lead to potential connector failure. Inspect the catheter frequently for nicks, scrapes and cuts, etc. which could impair catheter performance.
- The selection of the appropriate catheter length is at the sole discretion of the physician. To achieve proper tip placement, proper catheter length selection is important. Routine X-ray should always follow the initial insertion of the catheter to confirm proper placement prior to use.
- Early product failure might be as result of:
 - improper positioning of the catheter tip,
 - improper flushing and disconnection after use, or
 - aggressive insertion of male luer which may cause cracking of female luer.

Equipment required for venous access

- Sterile dressing
- Syringes and needles
- Facility for chest X-ray
- Sterile pack and antiseptic solution
- Appropriate central venous catheter
- Local anaesthetic (e.g. 5 ml lignocaine 1% solution)
- Saline or heparinised saline to prime and flush all lumens on the line after insertion
- Shaving equipment for the area if very hairy (especially if femoral access is required)
- Suture in case of fixation by suturing is determined (e.g. 2/0 silk on a straight needle)

GENERAL INSERTION TECHNIQUE FOR ALL ROUTES

- 1) Confirm that central venous access is needed and select the most appropriate route. Explain the procedure to the patient.
- 2) Shave the needle insertion area if very hairy.

- 3) Using a strict aseptic technique, prepare and check all the equipment for use. Ensure luer connector cap is removed from the distal (white) lumen to ensure the guidewire exits.
- 4) Prime the catheter lumens by flushing all of the lumens
- 5) Prepare the guidewire syringe by advancing the tip of the guidewire dispenser into the valve of the syringe. This mechanical action releases the silicone of the valve to receive the Nitinol guidewire
- 6) Prepare the skin and drape the area.
- 7) Infiltrate the skin and deeper tissues with local anaesthetic. In cases where difficulty is anticipated use the small local anaesthetic needle to locate the vein before using the larger needle. This reduces the risk of trauma to other structures.
- 8) Position the patient as for the specific route described – avoid long periods of head down, particularly in breathless patients.
- 9) Identify the anatomical landmarks for the chosen route and insert the needle at the recommended point. After the needle has penetrated the skin, aspirate gently whilst advancing the needle as directed until the vein is entered. Once the vein has been located, remove the syringe from the needle, if using conventional Seldinger technique, place thumb over needle to eliminate the risk of air embolism. If the vein is not found, slowly withdraw the needle whilst gently aspirating; often the vein has been collapsed and transfixed by the entry of the needle.
- 10) Advance the provided guidewire either through the guidewire syringe (modified Seldinger technique) or through the needle (Seldinger technique) into the vein. Monitor the marks on the guidewire to know the length of its advanced part. Generally, 20 cm for traditional Seldinger technique or 30 cm if using the guidewire introduction syringe.
- 11) If using the guidewire compatible syringe advance the guidewire dispenser into the syringe until it passes through the valve (Figure 1), then using the thumb on the wheel advance the guidewire all through the syringe to the inside of the attached introducer needle maintaining a firm grip on the syringe, while pushing the guidewire advancer tip into the syringe to an equivalent length to the desired position of the catheter tip.


Figure 1 – Guiding syringe

- 12) Once the guidewire has been advanced to the desired length to match the catheter tip, then proceed to the removal of the syringe. If using the guidewire syringe it is possible to remove needle, syringe and dispenser in one movement.
- 13) To prevent the guidewire migrating when the needle is being removed place pressure over the venopuncture whilst removing the needle.
- 14) Before inserting the dilator dip it into saline solution to activate the hydrophilic coating.
- 15) It may be necessary to dilate the hole in the skin. Use the scalpel, with the blade facing away from the guidewire, and make a small incision in the skin and fascia where the guidewire enters the patient. Thread the dilator over the wire into the vein with a twisting motion. Excessive force should not be needed. Remove the dilator taking care not to dislodge the guidewire. Ensure the dilator is not over-inserted (usually 4-6 cm maximum). Rotating the dilator may make insertion easier.
- 16) Thread the central venous catheter over the guidewire until the straight end of the wire protrudes out of LuerSafe® and whilst holding the wire still advance the catheter into the vein.
- 17) TAKE CARE not to allow the wire to be pushed further into the vein whilst advancing the catheter. Advance the catheter to the length determined anatomically or to the length determined through the use

of ECG lead if used. Remove the guidewire once the catheter is fully advanced. The catheter position should be confirmed with chest X-ray. Rotating the catheter may make insertion easier.

- 18) Aspirate blood up to the syringe to check that blood can be aspirated freely from all lumens of the catheter, flush with saline following push pause flushing technique and a positive disconnection technique
- 19) Secure the catheter in place with the suture and cover with a sterile dressing. Tape any redundant tubing carefully avoiding any kinking or loops which may snag and pull out the catheter.
- 20) Connect catheter lumens to a bag of intravenous fluid.

WARNING

- If any resistance is felt, then the needle should be pulled out with the wire still inside and the procedure repeated. This reduces the risk of entangling of the guidewire or its end being cut off by the needle tip.
- Over advancement of guidewire can result in serious injuries or arrhythmias.
- An ECG lead attached to the end of the guidewire can confirm its tip location at the right atrium through monitoring of an elevation of the P wave.
- The valve of the guidewire syringe should be opened by the guidewire advancer tip. Do not attempt passing the guidewire before opening the valve by the advancer tip otherwise the wire may be kinked or destroyed.
- Over advancement of the catheter can result in serious injury or arrhythmias.
- Use the marking on the catheter to determine the advanced length.
- Use 10 ml or greater syringe to flush the catheter to reduce risk of exceeding the pressure capacity of the catheter. If resistance is felt during flushing, no further attempt should be made. Further flushing could result in catheter rupture with possible leakage or immobilisation.

FIXATION BY SECONDARY FIXATION WING

- 1) Place the white part of the secondary fixation wing and place it as close as possible to the venopuncture.
- 2) Snap the blue part of the moveable wing over the white wing.
- 3) Suture the wings through the holes to the skin of the patient.


Figure 5

PRECAUTION

- Always ensure there is secure fixation of the catheter to the patient.
- Always ensure that the fixation is not causing catheter obstruction.

CHECKS BEFORE USING THE CATHETER

- 1) Ensure that fluid runs in freely and that blood draws back freely.
- 2) Take a chest X-ray (ideally erect) to check the position of the catheter tip and to exclude a pneumo-, hydro- or haemothorax. An early radiograph may not show up abnormalities and it may be best to wait 3-4 hours unless symptoms develop.

WARNING

- Always be sure that the catheter tip is in the superior vena cava and has not entered into the right atrium to avoid possible arrhythmias or mural injury.

PRACTICAL PROBLEMS COMMON TO MOST INSERTION TECHNIQUES

The table below lists some problems that may occur with any of the routes for central venous access.

Table 2 – Problems during central venous cannulation.

Problem type	Description
Arterial puncture	Usually obvious but may be missed in a patient who is hypoxic or hypotensive. Withdraw the needle and apply firm direct pressure to the site for at least 10 minutes or longer if there is continuing bleeding. If there is minimal swelling then retry or change to a different route.
Suspected pneumothorax	If air is easily aspirated into the syringe (note that this may also occur if the needle is not firmly attached to the syringe) or the patient starts to become breathless. Abandon the procedure at that site. Obtain a chest radiograph and insert an intercostal drain if confirmed. If access is absolutely necessary, then try another route ON THE SAME SIDE or either femoral vein. DO NOT attempt either the subclavian or jugular on the other side as bilateral pneumothoraces are produced.
Arrhythmias during the procedure	Usually from the catheter or wire being inserted too far (into the right ventricle). The average length of catheter needed for an adult internal jugular or subclavian approach is 15cm. Withdraw the wire or catheter if further than this.
Air embolus	This can occur, especially in the hypo-volaemic patient who is breathing spontaneously, if the needle is left in the vein whilst open to the air. It is easily prevented by ensuring that the patient is positioned in the Trendelenburg position (for jugular and subclavian routes) and that the guidewire is passed down the needle promptly.
The wire will not thread down the needle	Check that the needle is still in the vein. Flush it with saline. Try angling the needle so the end of it lies more along the plane of the vessel. Carefully rotate the needle in case the end lies against the vessel wall. Reattach the syringe and aspirate to check that you are still in the vein. If the wire has gone through the needle but will not pass down the vein it should be very gently pulled back. If any resistance is felt, then the needle should be pulled out with the wire still inside and the procedure repeated. This reduces the risk of the end of the wire being cut off by the needle tip.
Persistent bleeding at the entry	Apply firm direct pressure with a sterile dressing. Bleeding should usually stop unless there is a coagulation abnormality. Persistent severe bleeding may require surgical exploration if there is an arterial or venous tear.

POSSIBLE COMPLICATIONS

The table below lists some of the possible complications that can occur during the early and late stages of central venous cannulation.

Table 3 – Potential early and late complications

Early	Late
Arterial puncture	Air embolism Venous thrombosis
Bleeding	Catheter embolus Cardiac perforation and tamponade
Cardiac arrhythmias	Pneumothorax Infection

Early		Late
Injury to the thoracic duct	Injury to surrounding nerves	Hydrothorax

PREVENTION AND TREATMENT OF CATHETER DYSFUNCTION

Causes of early catheter dysfunction

- Mechanical compression (pinch off syndrome in subclavian catheter)
- Malposition of catheter tip
- Kinks
- Catheter migration
- Occlusion of side holes due to clotting or fibrin sheath formation or drug precipitation (some antibody locks or IV IgG)
- Patient position, especially if the catheter has not been well positioned and secured
- Loss of catheter integrity by infection

Methods that should be used to treat dysfunctional or non-functional catheter include

- Repositioning of a malpositioned catheter.
- Use of thrombolytics, as per hospital protocol.
- All catheter-related infections, except for catheter exit-site infections, should be addressed by initiating parenteral treatment with an antibiotic(s) appropriate for the organism(s) suspected. Definitive antibiotic therapy should be based on the organism(s) isolated.
- Catheters should be exchanged as soon as possible and within 72 hours of initiating antibiotic therapy in most instances, and such exchange does not require a negative blood culture result before the exchange. Follow-up cultures are needed 1 week after cessation of antibiotic therapy.

CARE OF CENTRAL VENOUS CATHETER

- Use an aseptic technique when inserting the catheter and any subsequent injections or changing fluid lines.
- When accessing LuerSafe wipe it in accordance with local hospital protocol but never use disinfectants with alcohol content 70% or higher and neither ethanol or propanol alcohols.
- Keep the entry site covered with a dry sterile dressing.
- Ensure the line is well secured to prevent movement (this can increase risks of infection and clot formation). It is also recommended that the tip position is monitored to ensure correct positioning.
- Change the catheter if there are signs of infection at the site.
- Remember to remove the catheter as soon as it is no longer needed. The longer the catheter is left in, the greater the risks of sepsis and thrombosis.
- It is recommended changing according to hospital protocols to reduce the risks of catheter related sepsis and thrombosis. However, providing that the catheter is kept clean (sterile injections and connections) and there are no signs of systemic sepsis, routine replacement may not be necessary. Repeated cannulation to change lines on a routine basis, rather than based on clinical need, can increase the risks to the patient.
- When changing the dressing, if using an alcohol based cleaning solution, ensure the site is dry by either waiting for 20 to 30 seconds or by drying with a clean, sterile swab before applying a dressing.
- For lumens not in use check at least once a day their patency by aspirating the lumen and flushing it with saline solution using a push-pause flushing followed by a positive disconnection.

CATHETER REMOVAL

Remove any dressing and suture material. It is recommended that the patient is in a supine or Trendelenburg position. Ask the patient to take a breath and fully exhale. Remove the catheter with a steady pull while the patient is holding their breath and apply firm pressure to the puncture site for at least 5 minutes to stop the bleeding. Excessive force should not be needed to remove the catheter. If it does not come out, try rotating it whilst pulling gently. If this still fails, cover it with a sterile dressing and ask an experienced person for advice.

CATHETER DISPOSAL

Used catheters should be disposed of in a sanitary container or according to hospital protocol to prevent possible contamination and cross infection.

DESCRIPTION OF THE MARKING SYSTEM

In order to achieve the required length, the catheter tube is marked with numerical numbers every 5 cm with dots every 1 cm in between, apart from the first 5 cm which is not marked.


5.....10.....15.....20.....25.....

Illustration purposes only, not to scale.

CODING CONVENTION

Altius Classic Central Venous Catheter Kits included in this range are coded as "KCS1XXYYZ", where:

- KCS1 stands for first generation uncoated;
- XX stands for the catheter length in centimetre;
- YY stands for the French size of the tube; and
- Z stands for the number of lumens of the catheter.

Advena Ltd.
Tower Business Centre, 2nd Floor,
Tower Street, Swatar, BKR 4013, Malta

Kimal PLC
Arundel Road, Uxbridge, Middlesex,
UB8 2SA, United Kingdom
Tel: +44 (0) 845 437 9540 Fax: +44 (0) 845 437 9541
www.kimal.co.uk

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