EN – SHORT TERM HAEMODIALYSIS CATHETER KIT	2
CZ – KATETRY PRO KRÁTKODOBOU HEMODIALÝZU	8
DA – KORTSIGTEDE HÆMODIALYSE KATETRE	14
DE – KURZZEIT-HÄMODIALYSE-KATHETER	20
EL – ΚΑΘΕΤΗΡΑΣ ΑΙΜΟΚΑΘΑΡΣΗΣ ΜΙΚΡΗΣ ΔΙΑΡΚΕΙΑΣ	26
ES – CATÉTERES PARA HEMODIÁLISIS TEMPORALES	32
FI – LYHYTAIKAINEN HEMODIALYYSIKATETRI	38
FR – CATHÉTERS COURT TERME POUR HÉMODIALYSE	44
IT – CATETERE PER EMODIALISI A BREVE TERMINE	50
NL – KORTE TERMIJN HEMODIALYSEKATHETERS	56
NO – KORTTIDS HEMODIALYSISK KATETRE	62
PL – CEWNIK KRÓTKOTERMINOWY DO HEMODIALIZY	68
PT – CATETERES DE CURTA DURAÇÃO PARA HEMODIÁLISE	74
SK – KATÉTRE NA KRÁTKODOBÚ HEMODIALÝZU	80
SL – HEMODIALIZNI KATETER ZA KRATKOTRAJNO UPORABO	86
SR – KATETER ZA KRATKOROČNU HEMODIJALIZU	92
SV – HEMODIALYSKATETRAR FÖR KORTVARIGT BRUK	98
TR – KISA SÜRELİ HEMODİYALİZ KATETERLERİ	104

Instructions for Use

INTENDED USE

Sterile single use device indicated for use in attaining short term (30 days) access for haemodialysis or aphaeresis.

DEVICE DESCRIPTION

The short term haemodialysis catheters are radiopaque and made from polyurethane catheter. The catheters are semi rigid but become softer at body temperature. Catheters tube profile can be single, dual or triple lumen. The catheter is tapered tipped with softer material. Catheter shape can be with both straight tube and extension, straight tube and J shape extension or pre-curved tube and straight extension. In "Smart type", annular groove at the side holes preventing their stick to the vessel wall.

How supplied

The device is sterilized by Ethylene Oxide. Contents are sterile and non-pyrogenic in unopened undamaged package. Do not use catheter if package has been damaged or has been opened.

Storage

Store at room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet lights. Rotate inventory so that catheters are used prior to expiration date on package label.

CONTRAINDICATIONS

The devices are contraindicated as follows:

- The catheter should not be placed in patient with bleeding disorders.
- When the presence of the other device related infection, bacteraemia or septicaemia is known or suspected.
- If severe chronic obstructive lung disease exists.
- Post irradiation of prospective insertion site.
- Previous episodes of venous thrombosis or vascular surgical procedure at the prospective placement site have occurred.
- Local tissue factors will prevent proper devices stabilization and/or access.

WHICH VEIN TO CANNULATE?

The preferred insertion site for dialysis catheters is the right internal jugular vein. Other options include the right external jugular vein, left internal and external jugular veins, subclavian veins and femoral veins. Subclavian access should be used only when no other upper-extremity or chest-wall options are available.

Table 1 – Patient evaluation prior to access placement

Consideration	Relevance
History of previous CVC	Previous placement of a CVC is associated with central venous stenosis.
Dominant arm	To minimize negative impact on quality of life, use of the non-dominant arm is preferred.
History of pacemaker use	There is a correlation between pacemaker use and central venous stenosis.

Consideration	Relevance
History of severe CHF	Placement of catheter may alter haemodynamic and cardiac output.
History of arterial or venous peripheral catheter	Previous placement of an arterial or venous peripheral catheter may have damaged target vasculature.
History of diabetes mellitus	Diabetes mellitus is associated with damage to vasculature necessary for internal accesses.
History of anticoagulant therapy or any coagulation disorder	Abnormal coagulation may cause clotting or problems with haemostasis of access site.
Presence of co-morbid conditions, such as malignancy or coronary artery disease, that limit patient's life expectancy	Morbidity associated with placement and maintenance of certain accesses may not justify their use in some patients.
History of vascular access	Previously failed vascular accesses will limit available sites for access; the cause of a previous failure may influence planned access if the cause is still present.
History of heart valve disease or prosthesis	Rate of infection associated with specific access types should be considered.
History of previous arm, neck or chest surgery/trauma	Vascular damage associated with previous surgery or trauma may limit viable access sites.

METHOD OF INSERTION

General preparation to obtain haemodialysis access

The basic preparation and equipment those are required for venous cannulation are the same regardless of the route or technique chosen. Clinicians who insert dialysis catheter should be taught the technique by an experienced colleague. If this is not possible then the access routes associated with the fewest complications are the femoral vein.

PRECAUTION

- **Ultrasound should be used in the placement of catheters.**
- **The position of the tip of any central catheter should be verified by a radiological means.**
- **Do not use absolute alcohol or acetone based product on the catheter. 2% chlorhexidine or Iodine based solution is recommended as antiseptic solution.**
- **It is not recommended to use ointments on catheters as it may cause its degradation.**
- **Over tightening of catheter luer may lead to its failure.**

Equipment required for venous access

- Sterile pack and antiseptic solution
- Local anaesthetic – e.g. 5ml lignocaine 1% solution
- Appropriate catheter for age/purpose
- Syringes and needles
- Saline or heparinised saline to prime and flush the line after insertion
- Suture material in case of fixation by suturing is determined – e.g. 2/0 silk on a straight needle
- Sterile dressing
- Shaving equipment for the area if very hairy (especially the femoral)

GENERAL INSERTION TECHNIQUE FOR ALL ROUTES

- 1) Confirm that central venous access is needed and select the most appropriate route. Explain the procedure to the patient.
- 2) Shave the needle insertion area if very hairy.
- 3) Using a strict aseptic technique, prepare and check all the equipment for use.
- 4) Prepare the skin and drape the area.

- 5) Infiltrate the skin and deeper tissues with local anaesthetic. In cases where difficulty is anticipated use the small local anaesthetic needle to locate the vein before using the larger needle. This reduces the risk of trauma to other structures.
- 6) Position the patient as for the specific route described – avoid long periods of head down, particularly in breathless patients.
- 7) Identify the anatomical landmarks for the chosen route and insert the needle at the recommended point. After the needle has penetrated the skin, aspirate gently whilst advancing the needle as directed until the vein is entered. If the vein is not found, slowly withdraw the needle whilst gently aspirating; often the vein has been collapsed and transfixed by the entry of the needle.
- 8) Advance a guidewire (Seldinger technique), into the vein, flexible J-shape end first, then remove the needle.
- 9) Advance the guide wire equivalent length to the desired position of the catheter tip.
- 10) It is necessary to dilate up the hole in the vein. Make a small incision in the skin and fascia where the wire enters the patient. Thread the dilator over the wire into the vein with a twisting motion. Excessive force should not be needed. Remove the dilator taking care not to dislodge the guidewire.
- 11) Thread the catheter over the guide wire until the end of the wire protrudes from the end of the catheter and whilst holding the wire still advance the catheter into the vein. TAKE CARE not to allow the wire to be pushed further into the vein whilst advancing the catheter.
- 12) Check that blood can be aspirated freely from all lumens of the catheter and flush with saline.
- 13) Secure the catheter in place with the suture and cover with a sterile dressing. Tape any redundant tubing carefully avoiding any kinking or loops which may snag and pull out the catheter.
- 14) Connect catheter to a bag of intravenous fluid or flush both lumen with appropriate anti thrombotic.

PRECAUTION

- **If any resistance is felt, then the needle should be pulled out with the wire still inside and the procedure repeated. This reduces the risk of entangling of the guidewire or its end being cut off by the needle tip.**
- **Over advancement of guidewire can result in serious injuries or arrhythmias.**

CHECKS BEFORE USING THE CATHETER

- 1) Ensure fluid runs in freely and that blood flows freely back.
- 2) If available, take a chest X-ray (ideally erect) to check the position of the catheter tip and to exclude a pneumo, hydro or haemothorax. An early radiograph may not show up abnormalities and it may be best to wait 3-4 hours unless symptoms develop. The tip of the catheter should lie in within the right atrium.
- 3) Ensure that the patient will be nursed and their access can be supervised. Give appropriate written instructions regarding how, and what it is to be used for, and who to contact if there is a problem between dialysis sessions.

PRACTICAL PROBLEMS COMMON TO MOST INSERTION TECHNIQUES

Table 2 – Problems during haemodialysis cannulation

Problem type	Description
Arterial puncture	Usually obvious but may be missed in a patient who is hypoxic or hypotensive. Withdraw the needle and apply firm direct pressure to the site for at least 10 minutes or longer if there is continuing bleeding. If there is minimal swelling then retry or change to a different route.
Suspected pneumothorax	If air is easily aspirated into the syringe (note that this may also occur if the needle is not firmly attached to the syringe) or the patient starts to become breathless, abandon the procedure at that site. Obtain a chest radiograph and insert an intercostal drain if pneumothorax is confirmed. If central access is absolutely necessary, then try another route ON THE SAME SIDE or either femoral vein. DO NOT attempt either the subclavian or jugular on the other side in order to avoid bilateral pneumothoraces.
Arrhythmias during the procedure	Usually from the catheter or wire being inserted too far (into the right ventricle). The average length of catheter needed for an adult internal jugular or subclavian approach is 15cm. withdraw the wire or catheter if further than this.

Problem type	Description
Air embolus	This can occur, especially in the hypovolemic patient, if the needle or cannula is left in the vein whilst open to the air. It is easily prevented by ensuring that the patient is positioned head down (for jugular and subclavian routes) and that the guidewire or catheter is passed down the needle promptly.
The wire will not thread down the needle	Check that the needle is still in the vein. Flush it with saline. Try angling the needle so the end of it lies more along the plane of the vessel. Carefully rotate the needle in case the end lies against the vessel wall. Reattach the syringe and aspirate to check that you are still in the vein. If the wire has gone through the needle but will not pass down the vein it should be very gently pulled back. If any resistance is felt then the needle should be pulled out with the wire still inside, and the procedure repeated. This reduces the risk of the end of the wire being cut off by the needle tip.
Persistent bleeding at the entry	Apply firm direct pressure with a sterile dressing. Bleeding should usually stop unless there is a coagulation abnormality. Persistent severe bleeding may require surgical exploration if there is an arterial or venous tear.

COMPLICATIONS

Table 3 – Potential early and late complications

Early		Late
Arterial puncture	Injury to surrounding nerves	Venous thrombosis
Bleeding	Air embolism	Cardiac perforation and tamponade
Cardiac arrhythmias	Catheter embolus	Infection
Injury to the thoracic duct	Pneumothorax	Hydrothorax

CONNECTION TO DIALYSIS MACHINE

Catheter is to be connected to the blood line of the dialysis machine. The blood line is a set of arterial and venous lines. After swabbing the female luer end, push straight with twisting the male luer of the blood tube inside the respectable catheter female luer end. Open the pinch clamp.

Ensure that the patient will be nursed during dialysis. Give appropriate written instructions regarding how, and what it is to be used for, and who to contact if there is a problem The catheter should allow a free flow of fluids. The free flow is usually indicated by flow of blood within the accepted venous and arterial pressure in the extracorporeal circuit of the dialysis machine.

DISCONNECTION FROM DIALYSIS MACHINE

Clamp the catheter extension line. Pull out with twisting the male luer of the blood tube from catheter female luer end. Connect the flushing device to the catheter female luer end then re- open the clamp. Flush the catheter according to hospital protocol. Clamp again the catheter and remove the flushing device. Use of dry gauze dressing combined with skin disinfection, using either chlorhexidine or povidone iodine solution, followed by povidone iodine ointment or mupirocin ointment at the catheter exit site are recommended at the end of each dialysis session.

PRECAUTION AGAINST MISUSE

Can be due to:

- Improper positioning of the catheter tip.
- Misconnection of catheter extension line(s) by connecting the venous line of the dialysis blood line to catheter extension line that is dedicated to the arterial line and marked red. This can lead to high recirculation rate and inefficient dialysis.
- Improper heparinization during the dialysis may result in blood clotting and obstruction of the catheter.
- Improper heparinization of the catheter between dialysis may result in thrombus formation.
- Inserting male luer aggressively may cracks catheter female luer.

PREVENTION AND TREATMENT OF CATHETER DYSFUNCTION

Catheters should be evaluated when they become dysfunctional. Dysfunction is defined as failure to attain and maintain an extracorporeal blood flow of 300 ml/min (for adult size catheter) or greater at a pre pump arterial pressure more negative than -250 mm Hg.

SIGNS OF CATHETER DYSFUNCTION – ASSESSMENT PHASE

- Blood pump flow rates < 300 ml/min
- Arterial pressure increases (< -250 mm Hg)
- Venous pressure increases (> 250 mm Hg)
- Conductance decreases (< 1.2): the ratio of blood pump flow to the absolute value of pre-pump pressure
- Unable to aspirate blood freely (late manifestation)
- Frequent pressure alarms – not responsive to patient repositioning or catheter flushing

PRECAUTION

- Trend analysis of changes in access flow is the best predictor of access patency and risk for thrombosis.

CAUSES OF EARLY CATHETER DYSFUNCTION

- Mechanical compression (pinch off syndrome in subclavian catheter)
- Malposition of catheter tip
- Kinks
- Catheter migration
- Side holes' occlusion due to clotting or fibrin sheath formation or stuck to vein wall
- Drug precipitation (some antibody locks or IV IgG)
- Patient position especially in not well fixed and secured catheter
- Loss of catheter integrity by infection

METHODS THAT SHOULD BE USED TO TREAT DYSFUNCTIONAL OR NON-FUNCTIONAL CATHETER INCLUDE

- Repositioning of a malpositioned catheter.
- Change patient position, ask him to cough or vigorous flush (if no resistance is felt) trying to dislodge side holes a way from vein wall.
- Fibrin sheath stripping using a snare if a fibrin sheath is present.
- Exchanging the thrombosed catheter over a guidewire, if a fibrin sheath is present or if the catheter is malpositioned or of inadequate length.
- Use of thrombolytics, as per hospital protocol.
- Treatment of an infected HD catheter should be based on the type and extent of infection.
- All catheter-related infections, except for catheter exit-site infections, should be addressed by initiating parenteral treatment with an antibiotic(s) appropriate for the organism(s) suspected.
- Definitive antibiotic therapy should be based on the organism(s) isolated.
- Catheters should be exchanged as soon as possible and within 72 hours of initiating antibiotic therapy in most instances, and such exchange does not require a negative blood culture result before the exchange. Follow-up cultures are needed 1 week after cessation of antibiotic therapy.

At the end of the dialysis session, clamp catheter extension line(s) and cap the catheter with the injection cap. This should be followed by injecting heparin or other anti-thrombotic at each lumen (according to lumen priming volume) in the catheter via the injection caps.

CARE OF CATHETER BETWEEN DIALYSIS

- Insertion site should be inspected for possible bleeding.
- Anti-thrombotic should be regularly injected to the catheter to prevent catheter thrombus and obstruction.

CATHETER REMOVAL

Remove any dressing and suture material. Ask the patient to take a breath and fully exhale. Remove the catheter with a steady pull while the patient is breath holding and apply firm pressure to the puncture site for at least 5 minutes to stop the bleeding. Excessive force should not be needed to remove the catheter. If it does not come out, try rotating it whilst pulling gently. If this still fails, cover it with a sterile dressing and ask an experienced person for advice.

CATHETER DISPOSAL

Used catheter should be disposed in sanitary container to prevent possible contamination and cross infection.

DESCRIPTION OF MARKING SYSTEM

The catheter tube is marked for effective length in numerical number every 5 centimetres and dot every one centimetre however the first 5 cm is not marked.

5.....10.....15.....20.....25.....

Illustration purposes only, not to scale.

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Advena Ltd.